



Topcon RM-A7000 ARK

Reconditioned • Calibrated • Patient Ready



This compactly designed unit features a built-in TV monitor to facilitate alignment and measurement: an automatic PD measurement: an automatic sleep mode and a built-in printer.

Up to 10 reading per eye can be taken. With three or more readings and average value is calculated and displayed on the print-out. Automatic right and left eye determination. Fixation target incorporates an auto-fog system. RS-232 interface allows for online transfer of measurement of data.

For a Bigger Return on Your Investment

Call, text, WhatsApp: 949-362-1150

info@usedeyeequipment.com