

PRE-OWNED EXCIMER/FEMTO LASERS

Refurbished, like new, patient ready!



Allegretto Wavelight Excimer Laser 400 Hz

WaveFront Optimization: Customized treatment based on the unique characteristics of the eye.



AMO Visx Star S4 Excimer Laser

Wavescan and Iris Registration (IR). Some S3 and solo units are available.



Wieland EC 5000 Excimer Laser Board AGS units are available options, eye tracker, custom



TruLase FEMTO 50 and FS Laser claims to delivers micron-level accuracy 100% greater^{1,2}



Ziemer Femtosecond Laser performed with the Ziemer FEMTO LDV Femtosecond Surgical
Contact eye equipment 0011493621450