

Topcon NW-200 Non Mydriatic Fundus Camera

Written by admin



Topcon NW-200 Non Mydriatic Fundus Camera

Written by admin

Topcon NW-200 Non-Mydriatic Fundus Camera

Refurbished, Calibrated, Reliable, Guaranteed

The TRC-NW200 Telemedicine System is simple and easy to use, making it an ideal tool for the health professional looking for a quick and accurate way to screen and document retinal pathology.

It combines ImageNet™ Lite, the recognized standard in color fundus imaging software, with the Topcon TRC-NW200 Non-Mydriatic Retinal Camera. Combined with a laptop PC, the TRC-NW200 System affords a portable screening solution. It can also be mated with a Tower PC.

For a Bigger Return on Your Investment

Call, text, WhatsApp: 949-362-1150

info@usedeyeequipment.com

Topcon NW-200 Non Mydriatic Fundus Camera

Written by admin
