

PRE-OWNED FUNDUS CAMERAS

Refurbished, like new, patient ready!



Canon CR1/ Digital Non-Mydriatic Camera

The CR-1 Mark II can achieve extremely detailed, high-resolution diagnostic images of the retina and enables a full 45-degree view angle to help the eye care professional detect and monitor ocular conditions. When attached to a Canon EOS 50D Digital SLR Camera which provides up to 15.1 million pixel resolution, the CR-1 Mark II is capable of delivering clear and detailed diagnostic images for virtually immediate review.



Canon CR-2 Digital Retinal Camera. Setup and alignment conditions is easier



Topcon TRC-NW900 Digital Retinal Camera. Topcon TRC-NW900 is ideal for patients, and can be used for fundus photography.



Zeiss Visucam 700 Digital Retinal Camera. Eye Disease



Full-Frame Digital Retinal Camera with High-Speed Image Transfer and Storage, with
info@usedeyequipment.com